


# Multi Purpose Safety Station



The Multi Purpose Safety Station features:

1. Eye/Face Wash
2. Drench Shower/ Face Wash
3. Safety Storage Cabinet



- All in only 6 sq. ft. of Floor Space!
- Constructed of 2" panels for durability, sanitation & corrosion resistance.
- Shipped completely assembled, ready to install.
- Fixtures meet ANSI Z358.1 standard for Safety Showers.
- Ideal for work sites where persons may be exposed to potentially hazardous chemicals.
- The HEMCO Safety Station is equipped with a Drench shower, and deck mounted eye face wash.
- Towels and first aid supplies can be stored in the cabinet below.
- Safety Stations are fully assembled and ready for installation to water supply and drain.
- Supply 1.5" NPT, Drain 1.25" NPT.

## Safety Station Specifications:

Unit is equipped with Drench shower, and deck mounted eye/face wash & storage cabinet for safety supplies.

30" wide Safety Station

Cat. No. 78001

36" wide Safety Station

Cat. No. 78002

