

# C - Frame Table

The HEMCO C - Frame Table is designed to support Demonstration Hoods, VSE Enclosures and light weight Laboratory Equipment.



- **Mobile Table:** Sizes from 30", 36" & 48" wide
  - 36" Standard Height
  - Constructed of 1" X 2", 16 gauge rectangular tubing, with welded construction
  - (4) adjustable leg levelers included
  - Epoxy Resin Work Surface measures 3/4" thick X 24" deep X Length depending upon size
- Cat. No.**
- Powder Coated finish resists chemicals  
Color: White

## Specifications:

HEMCO C-Frame table is designed to hold:

- MicroFlow I Workstations
- MicroFlow II Workstations
- Demonstration Hoods
- VSE Vented Safety Enclosures
- Table Top Fume Hoods

Cat. No.	Width A	Height B	Depth C
50203	30"	36"	24"
50303	36"	36"	24"
50403	48"	36"	24"

