

Click to Clear Form

Clear

caution: this clears all form fields

HEMCO Corporation ISO 9001: 2008

Click to Submit Form

Submit

window opens, select your e-mail service, add a message, press send

UniMax Specification Planner

HEMCO Corporation PH:816-796-2900 Fax:816-796-3333 info@HEMCOCORP.COM

UniMax Specification Planner

CATALOG#

SHEET 1 OF 3

CUSTOMER BY PHONE FAX EMAIL DATE

CUSTOMER APPROVAL APPROVED AS SHOWN APPROVED AS NOTED NOT APPROVED

CUSTOMER SIGNATURE DATE




Table with columns: Dimensions (Width A, Height B), Depth C (48" outside, 72" outside), and Number of Exhaust Outlets. It lists various size combinations and their corresponding catalog numbers.

Select a UniMax hood size. A. select height, B. select width, C. select depth

A) ROOM DESIGN PARAMETERS (attach floor plan if possible)

- 1. Hood Size: (A) Width X (B) Depth X (C) Height, Catalog No.
2. Internal Dimension requirements (A) Width X (B) Depth X (C) Height
3. Process to be performed and type of equipment in hood
4. Airflow requirements in room:

B) DOORS

- 1. Doors: (FG) Framed Glass with overhead suspension (UG) Unframed Glass with overhead suspension (FS) Framed Glass with Floor Sliding Track (SC) Strip Curtains
2. Access: Front Rear Right Left

Click to Clear Form

Clear

caution: this clears all form fields

HEMCO Corporation ISO 9001: 2008

Click to Submit Form

Submit

window opens, select your e-mail service, add a message, press send

UniMax Specification Planner

HEMCO Corporation PH:816-796-2900 Fax:816-796-3333 info@HEMCOCORP.COM

UniMax Specification Planner

CATALOG#

SHEET 2 OF 3

CUSTOMER \_\_\_\_\_
BY \_\_\_\_\_ DATE \_\_\_\_\_
PHONE \_\_\_\_\_
FAX \_\_\_\_\_
EMAIL \_\_\_\_\_

CUSTOMER APPROVAL [ ] APPROVED AS SHOWN
[ ] APPROVED AS NOTED (PLEASE MARK CORRECTIONS IN RED)
[ ] NOT APPROVED (PLEASE INDICATE REASON)

CUSTOMER SIGNATURE \_\_\_\_\_ DATE \_\_\_\_\_

C) LIGHTING / ELECTRICAL

1. Electrical service in building: 440V / 480V. 30 \_\_\_\_\_ 220V /240V.30 \_\_\_\_\_

110V / 208V.30 \_\_\_\_\_ is 227V available for lighting? Yes [ ] No [ ] Not Sure [ ]

2. Lighting: Fluorescent \_\_\_\_\_ UV \_\_\_\_\_

Explosion Proof Incandesant \_\_\_\_\_ None \_\_\_\_\_

3. Electrical plug in requirements \_\_\_\_\_ Standard \_\_\_\_\_ Explosion Proof \_\_\_\_\_

4. Cord Port \_\_\_\_\_ Location \_\_\_\_\_

D) PLUMBING SERVICES

1. D.I Faucet (PVC) \_\_\_\_\_ D.I. Faucet (chrome) \_\_\_\_\_ D.I. Faucet (tin lined) \_\_\_\_\_

Water Faucet (cold water) \_\_\_\_\_ Service fixture (turret) \_\_\_\_\_ Remote Service fixture \_\_\_\_\_

(specify Air, Gas, Water, Vaccuum etc.) \_\_\_\_\_

2. Wall Mount Oval Cup Sink \_\_\_\_\_

3. Fume Hood Fire Extinguisher \_\_\_\_\_

E) EXHAUST / FILTRATION

1. Exhaust Collar Diameter \_\_\_\_\_

2. Exhaust Fan: Yes \_\_\_\_\_ No \_\_\_\_\_ If Yes, CFM \_\_\_\_\_ Or ACH \_\_\_\_\_

3. Length of Exhaust Run \_\_\_\_\_ Number of elbows \_\_\_\_\_

F) WORK SURFACES

1. Epoxy Resin (ER-WS) \_\_\_\_\_ Phenoloc Resin (PR-WS) \_\_\_\_\_ Stainless Steel (SS-304) \_\_\_\_\_

G) DISTILLATION APPARATUS GRID

1. Aluminum Rod assembly \_\_\_\_\_ Width \_\_\_\_\_

2. Stainless Steel Rod assembly \_\_\_\_\_ Width \_\_\_\_\_

3. Plastic (FG) Rod assembly \_\_\_\_\_ Width \_\_\_\_\_

Click to Clear Form

Clear

caution: this clears all form fields

HEMCO Corporation ISO 9001: 2008

UniMax Specification Planner

Click to Submit Form

Submit

window opens, select your e-mail service, add a message, press send

HEMCO Corporation PH:816-796-2900 Fax:816-796-3333 info@HEMCOCORP.COM

UniMax Specification Planner


CATALOG#

SHEET 3 OF 3

CUSTOMER \_\_\_\_\_
BY \_\_\_\_\_ DATE \_\_\_\_\_
PHONE \_\_\_\_\_
FAX \_\_\_\_\_
EMAIL \_\_\_\_\_

CUSTOMER APPROVAL [ ] APPROVED AS SHOWN
[ ] APPROVED AS NOTED (PLEASE MARK CORRECTIONS IN RED)
[ ] NOT APPROVED (PLEASE INDICATE REASON)

CUSTOMER SIGNATURE \_\_\_\_\_ DATE \_\_\_\_\_



For Containment basin add suffix CB to UNIMAX model #

Table with 3 columns: Dimensions (Width A), Depth C (48" and 72"), and model numbers for various basin sizes (6', 8', 10', 12', 16', 20', 24').

H) HAZMAX SPILL CONTAINMENT BASIN

1. Basin Size: (A) Height \_\_\_\_\_ X (B) Width \_\_\_\_\_ X (C) Depth \_\_\_\_\_, Catalog No. \_\_\_\_\_

2 Volume of liquid to be contained \_\_\_\_\_